



# BE SAFE

## PARAMEDICAL

WWW.BE-SAFE.CO.ZA | INFO@BE-SAFE.CO.ZA

# Cri+ Care<sup>®</sup>

"Progress Saves Lives..."



## Burn-Eaz<sup>®</sup> Mini Grabber Bag

Product Code: REF:BAG.GM.BE.PVC.N

The baby brother of the Grabber Bag, for those requiring slightly less packing volume than the Grabber Bag.

- Manufactured from waterproof 550G PVC.
- YKK zips.
- 5 internal removable storage pouches.
- Storage pouches have transparent covers for quick content identification made from 250 micron PVC.
- Double handles make for comfortable carrying.

Weight: 0.58kg

Dimensions: 32cm x 25cm x 9cm



An ISO 9001:2015 & 13485:2016 Certified Company  
[www.be-safe.co.za](http://www.be-safe.co.za)  
[info@be-safe.co.za](mailto:info@be-safe.co.za)

CAPE TOWN: +27 21 788 4681  
JOHANNESBURG: +27 11 794 7476  
DURBAN: +27 31 702 3219