



BE SAFE

PARAMEDICAL

CritiCare[®]

"Progress Saves Lives..."



WWW.BE-SAFE.CO.ZA | INFO@BE-SAFE.CO.ZA

CritiCare[®] AmpPAC[™] Pro Micro Drug Bag

Product Code: REF:BAG.AMP.PRO.M.D



The newest and tiniest member of our CritiCare[®] Drug Bag Family, the AmpPAC[™] Pro may be just the perfect drug bag for you.

- Manufactured from high tenacity laminated, water and UV resistant 420D nylon with a PVC backing.
- YKK zips.
- Central divider for storage of syringes and needles.
- Carry Handles.
- 10 small foam mesh removable pouches with clear window for easy contents identification.
- One medium size removable pouch with Clear window for easy contents identification.

Weight: 0.33kg

Dimensions: 24cm x 17cm x 5.5cm

An ISO 9001:2015 & 13485:2016 Certified Company
www.be-safe.co.za
info@be-safe.co.za

CAPE TOWN: +27 21 788 4681
JOHANNESBURG: +27 11 794 7476
DURBAN: +27 31 702 3219