



BE SAFE

PARAMEDICAL

WWW.BE-SAFE.CO.ZA | INFO@BE-SAFE.CO.ZA

Cri+ Care[®]

"Progress Saves Lives..."



CritiPack™ Bougie-aided Surgical Cricothyroidotomy

Product Code: REF:CRITI.P.BSC

CritiPack[®] Bougie-aided Surgical Cricothyroidotomy Kit provides access to all the equipment required to perform rapid surgical cricothyroidotomy in the field

- 1 x Pair XL Black Vitriile Gloves
- 1 x Disposable Scalpel
- 1 x 6.0mm Cuffed Tracheostomy Tube
- 1 x 10ml syringe with Luer Lock
- 1 x 15Fr Intubating Bougie
- 4 x Alcohol Swabs
- 1 x Pack of 5 Sterile 100 X 100 X 8PLY Gauze Swabs
- 1 x 400 x 400 protection sheet
- Kit packaged in sterile pouch

Weight: 0.2 kg

Dimensions: 34cm x 14cm x 4cm