



BE SAFE

PARAMEDICAL

WWW.BE-SAFE.CO.ZA | INFO@BE-SAFE.CO.ZA

CritiCare®

"Progress Saves Lives..."



CritiCare® Armsling

Product Code: REF:AS.CRITI.BLK

The CritiCare® Arm Sling is designed to immobilize injured shoulders and elbows during the rehabilitation phase of care

- Manufactured from high quality foam mesh for comfort and improved airflow
- Universal and fully adjustable shoulder sling
- Elastic Thumb Loops
- Padded shoulder strap
- Fits most sizes
- Hand washable with mild soap in warm water. Rinse well and air dry. Do not bleach, dry clean or tumble dry

Weight: 0.14 kg

Dimensions: 42cm x 18cm x 2cm



An ISO 9001:2015 & 13485:2016 Certified Company
www.be-safe.co.za
info@be-safe.co.za

CAPE TOWN: +27 21 788 4681
JOHANNESBURG: +27 11 794 7476
DURBAN: +27 31 702 3219