

Reusable Body Bags

VARIOUS TYPES

Cri+Care[®]

"Progress Saves Lives..."



Reusable body bags are used for storing a body and makes it easier to move a body from one place to another. It's also easy to store cause it can be folded up to take less space.

Variations, Product Codes & Item Dimensions:

Product Info	Order Code	Dimensions Folded	Dimensions Unfolded	Weight
Adult	REF:BAG.BODY.A	40 x 31 x 9 cm	220 x 94 x 1 cm	1 kg
Child	REF:BAG.BODY.C	31 x 16 x 4 cm	122 x 66 x 1 cm	275 g

Features:

- Comes with handles
- Includes Zips
- Made from heavy duty material
- Leak proof and Waterproof

Be Safe Paramedical cc South Africa

Cri+Care[®]
"Progress Saves Lives..."

An ISO 9001:2015 & 13485:2016 Certified Company
www.be-safe.co.za
info@be-safe.co.za

CAPE TOWN: +27 21 788 4681
JOHANNESBURG: +27 11 794 7476
DURBAN: +27 31 702 3219

V01-2020